

High and Low Demand Analysis

Dreamy Joe

Date: 3/8/2004 **Time:** 9:35:00 AM **Age:** 9 years 7 months **On Meds:** N
Diagnosis: **Group Code:** **ID Code:**
Medications:
Comment:
End Note:

Auditory Response Validity Check = Valid
Visual Response Validity Check = Valid

	Quotient Scores		
	Combined	Auditory	Visual
HIGH DEMAND TEST CONDITIONS			
-----	-----	-----	-----
COMPETENCE (Impaired-Superior)	43	45	52
-----	-----	-----	-----
Prudence (Impulsive-Prudent)	98	96	99
Quickness (Slow-Quick)	80	81	83
Steadiness (Faltering-Steady)	36	63	34
Stability (Erratic-Stable)	52	51	70
LOW DEMAND TEST CONDITIONS			
-----	-----	-----	-----
MAINTAINABILITY (Impaired-Superior)	33	57	24
-----	-----	-----	-----
Acuity (Unobservant-Sharp)	39	71	28
Dependability (Inconstant-Resolute)	58	69	62
Reliability (Careless-Reliable)	95	99	91
Swiftiness (Slack-Swift)	70	77	74

The global scales are capitalized and are based on the combined quotient scores of their respective individual scales. Thus, no raw data exists for them.

Note: Quotient scores which are between zero and 59 reflect extreme deficits.